



STATCLEAN TECHNOLOGY (S) PTE. LTD.

30 Kaki Bukit Road 3, # 03-12/13/14, Empire Technocentre, Singapore 417819.

Tel: 65-6746 2706 Fax: 65-6746 3045 / 65-6746 3755

Email: sales@statclean.com Website: www.statclean.com

Sticky Mat, Grey

(STMA-SK)

Features

- 99.8 % of dust collection at several steps on the mats.
- Collected dust on the mats is not migrated.
- Comfortable walking and smooth passing of hand carts.
- Antimicrobial sticky mat is available for pharmaceutical, biological, medical areas.

Size: 18" x 36", 18" x 45", 24" x 36", 26" x 45", etc.

Colour: Grey

Packing: 30 Sheets/Mat, 8 Mats / Carton

Thickness: 0.045 mm

Tacky Level : 190/250/300 ± 50 (Top)

(gf/25mm) 250/120/120 ± 100(Back)

700 ± 100(Bottom)

	Lists of fungi to be used	Result	Test method
D, E Antimicrobial Test	Aspergillus niger(ATCC 9642)	No growing	ASTM G 21
	Chaetomium globosum(ATCC 6205)		
	Penicillium pinophilum(ATCC 11797)		
	Gliocladium virens(ATCC 9645)		
	Aureobasidium pullulans(ATC C 15233)		